

Enzyme activity (10 families)

Supplies and Reagents:

NaHPO₄ buffer: pH 7, 100 mM

Substrate: Ask Jie

Cofactor: prepare the cofactor (Table 1)

Table 1. Amount of buffer, substrate, cofactors and H₂O in each reaction

enzyme family	pH buffer	substrate (final conce	cofactors (final concentration)			H ₂ O
PAL						
	100mM	50uM				add 25 ul to 100ul
	50 ul	5 ul				
4CL			CoA	ATP	MgCl ₂	
	100mM	50uM	100uM	2 mM	2mM	add 16 ul to 100ul
	50 ul	5 ul	5 ul	2 ul	2 ul	
C4H, C3H, Cald5H			NADPH	G6P	G6PdeH	
	100mM	50uM	1mM	10mM	2ul	add 4 ul to 100ul
	50 ul	5 ul	4 ul	10 ul	2 ul	
ComT, CoAOMT			MgCl ₂	SAM		
	100mM	50uM	2mM	250uM		add 20 ul to 100ul
	50 ul	5 ul	2 ul	3 ul		
CAD			NADPH	G6P	G6PdeH	
	100mM	50uM	1mM	10mM	2ul	add 20 ul to 100ul
	50 ul	5 ul	4 ul	10 ul	2 ul	
CCR			NADPH	G6P	G6PdeH	
	100mM	50uM	1mM	10mM	2ul	add 4 ul to 100ul
	50 ul	5 ul	4 ul	10 ul	2 ul	
HCT			shikimic acid			
	100mM	50uM	500 uM			add 24 ul to 100ul
	50 ul	5 ul	1 ul			

Important Notes:

Before HPLC, use the centrifuge for all the tubes in 20,000Xg for 30 min.

Protocol:

1. Enzyme reaction

1. For each enzyme family, add the buffer, substrate, cofactors and H₂O to 6 tubes (Table 1).
3 Tubes will used as control, and 3 tubes will used as reaction.
2. 3 Tubes will used as background, and add the buffer (50 ul) and H₂O (30 ul) to each tube.
3. 3 Tubes will used as reference for #22, and add the same cofactors as CAD.
4. 3 Tubes will used as reference for #24, and add the same cofactors as CAD.
5. Incubate all the tubes in 30 C.

6. Start the reaction by adding 20 ul xylem protein to reaction tubes. For control, background, and reference tubes, add 20 ul boiled xylem protein.
7. Stop the incubation of all tubes (except the CAD reaction tubes and reference tubes) at 30 min by adding 3M TCA (5 ul).
8. Stop the incubation of CAD and reference tubes at 30 min by heating the tubes in 99 C for 5 min.
9. Use the centrifuge for all the tubes in 20,000Xg for 30 min.
10. Transfer the supernatant to HPLC insert.

2. Preparation of the equipment

1. Check the column, in this analysis C18 column will be used.
2. Turn on the five sections of the equipment, click on Instrument online to set up for the analysis.
3. Check the level of the bottle before starting any sequence. Fill if necessary and reset the volume of the bottles.
4. A standard method for the enzyme activity named **Enzyme activity 1 on the HPLC machine** should be used.
5. Set up sequence and run.

3. Data analysis

1. Get the peak area by using Rt and wave length in Table 2.
2. Use the template to get the relative activity and specific activity of each enzyme.

Table 2. Rt and Wave length for each compound

enzyme family	substrate	Rt(min)	wave length (nm)
PAL	Phenylalanine		
C4H	Cinnamic acid	9.442	280
4CL	p-Coumaric acid	5.352	310
4CL	Caffeic acid	3.817	325
HCT	p-Coumaroyl-CoA	4.503	325
C3H	p-Coumaroyl shikimic acid	4.903	310
CCoAOMT	Caffeoyl-CoA	3.81	325
CCR	Feruloyl-CoA	4.709	325
CAld5H	Coniferaldehyde	8.385	325
COMT	5-Hydroxyconiferaldehyde	6.427	325
CAD/SAD	Coniferaldehyde	8.385	325
CAD/SAD	Sinapaldehyde	8.31	325
	p-Couamaryl alcohol	5.956	260
	Coniferyl alcohol	6.548	260
	Sinapyl alcohol	6.437	280

