

# Protocol for *P. trichocarpa* transformation

20120206 by Chenmin

## **Agrobacterium strain:**

*Agrobacterium tumefaciens* C58 w/ pBI121 construct

## **Plant material:**

5-6 month old *Populus trichocarpa*, 150-200 cm in height, stem segments (~25cm) between 5-9 internodes (depend on trees) from greenhouse

## **Transformation Protocol:**

1. Streak *Agrobacterium* single colony on a LB plate containing appropriate antibiotics. Incubate at 28°C for 2 days.
2. Pick up one isolated *Agrobacterium* colony to inoculate 3mL of LB medium with the appropriate antibiotics and grow for 20 hours in a shaker at 28°C, 200rpm.
3. *Agrobacterium* cells were grown for 20 hours in LB medium supplemented with the appropriate antibiotics on a shaker at 28 °C and 200 rpm. (Acetosyringone was added to LB medium )
4. Razor blade soaked in 70% ethanol for ~30 min.
5. Pick up trees 5-10 internodes in greenhouse and wash plant materials collected in 10 % Clorox solution for 20 minutes and at least three times rinse in sterile dd water.  
**Notes:** Keep in mind, pick up internodes varies depending on trees condition. After rinse immediately follow next step.
6. Stem segments (4 mm in length) were excised using razor blade and inoculated by swirling in an *Agrobacterium* culture for five minutes.
7. The inoculated explants were co-cultivated on CIM1 at 25 °C in darkness for two days.
8. Explants were washed four times in sterile deionized water, blotted dry on a sterile filter paper, and cultured in the dark on CIM2 at 25 °C. First, explants were subcultured twice every five days. Then, they were subcultured every 14 days until 45-60 days of co-cultivation.  
**Notes:** White calli should appear about 1 month in this stage, if not, prolong the time in the dark.

*Agrobacterium* control is crucial in this stage. Check your plate regularly. If infested, cleanup the infested explants, and transfer those good ones to a new plate.

9. Explants with putative transgenic callus were transferred in light at 25 °C onto SIM1 and subcultured every 14 days until 42 days.
10. Then explants with putative transgenic callus were inoculated on SIM3 for 1 ~ 4 months in light at 25 °C, subcultured every 20 days.
11. Putative transgenic shoots (1 cm in length) were cut off and inoculated on RM at 25 °C in low light. Whole plantlets with good root system were gained after approximately 1-1.5 months.
12. When rooted shoot reaches Magenta box top, you might confirmation by PCR and propagate copies. Later, you can move to green house with control wide type.

#### **Agrobacterium medium:**

LB, add kanamycin 50mg/L, genmycin 50mg/L, AS 20um

#### **Culture Medium**

##### CIM1

Basal medium: 1X MS salts and 1X vitamins

Kinetin: 0.5 mg/L

2,4-D: 0.05 mg/L

Inositol: 100mg/L

Sucrose: 30g/L

Adjust pH to 5.7 by 1N KOH, add

Agar: 6.5 g/L

Autoclaved, then add

Coconut water: 10 %

##### CIM2

Basal medium: 1X MS salts and 1X vitamins

Kinetin: 0.5 mg/L

2,4-D: 0.05 mg/L

Inositol: 100mg/L

Sucrose: 30g/L

Adjust pH to 5.7 by 1N KOH, add

Agar: 6.5 g/L

Autoclaved, then add filter sterilized

Kanamycin 50 mg/L

Cefotaxime 500 mg/L

##### SIM1

Basal medium: 1/2X MS salts and 1X vitamins

Inositol: 100mg/L  
Sucrose: 30g/L  
Adjust pH to 5.7 by 1N KOH, add  
Agar: 6.5 g/L  
Autoclaved, then add filter sterilized  
Thidiazuron: 0.132 mg/L  
Kanamycin 100 mg/L  
Cefotaxime 500 mg/L

### SIM3

Basal medium: 1X WPM  
Inositol: 100mg/L  
Sucrose: 30g/L  
Adjust pH to 5.7 by 1N KOH, add  
Agar: 6.5 g/L  
Autoclaved, then add filter sterilized  
Zeatin: 2.0 mg/L  
Kanamycin 50 mg/L  
Cefotaxime 500 mg/L

### RM

Basal medium: 1 X WPM  
Inositol: 100mg/L  
Sucrose: 20g/L  
Ca-gluconate: 650mg/L  
MES: 500mg/L  
Adjust pH to 5.75 by 1N KOH, add either  
Agar: 6.5g/L or phytogel: 3g/L  
Autoclaved